



ETF A rabbit pAb

□□	YP-Ab-18083
□□□	IgG
□□	WB
□□	Human;Mouse;Rat
□□	ETF A
□□□□	ETF A
□□□	Synthesized peptide derived from human ETF A
□□□	This antibody detects endogenous levels of ETF A at Human, Mouse,Rat
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	30 kDa
□□	Heterodimeric electron transfer flavoprotein that accepts electrons from several mitochondrial dehydrogenases, including acyl-CoA dehydrogenases, glutaryl-CoA and sarcosine dehydrogenase . It transfers the electrons to the main mitochondrial respiratory chain via ETF-ubiquinone oxidoreductase (ETF dehydrogenase) . Required for normal mitochondrial fatty acid oxidation and normal amino acid metabolism .
□□□□	Mitochondrion matrix .
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images