



TDO2 rabbit pAb

□□	YP-Ab-18119
□□□	IgG
□□	WB
□□	Human;Mouse;Rat
□□	TDO2
□□□□	TDO2 TDO
□□□	Synthesized peptide derived from human TDO2
□□□	This antibody detects endogenous levels of TDO2 at Human, Mouse,Rat
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	406 aa, 48 kDa
□□	Heme-dependent dioxygenase that catalyzes the oxidative cleavage of the L-tryptophan (L-Trp) pyrrole ring and converts L-tryptophan to N-formyl-L-kynurenine. Catalyzes the oxidative cleavage of the indole moiety.
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images