



CLMP rabbit pAb

□□	YP-Ab-18276
□□□	IgG
□□	WB
□□	Human;Mouse;Rat
□□	CLMP
□□□□	CLMP ACAM ASAM UNQ318/PRO363
□□□	Synthesized peptide derived from human CLMP □□□□
□□□	This antibody detects endogenous levels of CLMP at Human, Mouse,Rat
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	373 aa, 41 kDa
□□	May be involved in the cell-cell adhesion. May play a role in adipocyte differentiation and development of obesity. Is required for normal small intestine development.
□□□□	Cell junction, tight junction . Cell membrane ; Single-pass type I membrane protein .
□□□□	Predominantly expressed in epithelial cells within different tissues and in the white adipose tissue. Expressed at high levels in small intestine and placenta, at intermediate levels in the heart, skeletal muscle, colon, spleen, kidney and lung and at low levels in the liver and peripheral blood leukocytes. Highly abundant in the intestine during embryo and fetal development (at protein level).
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

