



SUN1 Rabbit pAb

Catalog No	YP-Ab-18412
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB
Gene Name	SUN1 KIAA0810 UNC84A
Immunogen	Synthesized peptide derived from human SUN1
Specificity	This antibody detects endogenous levels of SUN1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0810; Protein unc 84 homolog A; Sad1/unc 84 protein like 1; SUN1; unc 84 homolog A (C. elegans); UNC84A; SUN1 KIAA0810 UNC84A; SUN domain-containing protein 1 (Protein unc-84 homolog A) (Sad1/unc-84 protein-like 1)
Observed Band	90-100 kDa
Calculated Molecular Weight	90 kDa
Cell Pathway	Nucleus inner membrane ; Single-pass type II membrane protein . At oocyte MI stage localized around the spindle, at MII stage localized to the spindle poles. .
Tissue Specificity	
Function	As a component of the LINC (Linker of Nucleoskeleton and Cytoskeleton) complex involved in the connection between the nuclear lamina and the cytoskeleton . The nucleocytoplasmic interactions established by the LINC complex play an important role in the transmission of mechanical forces across the nuclear envelope and in nuclear movement and positioning (By similarity). Required for interkinetic nuclear migration (INM) and essential for nucleokinesis and centrosome-nucleus coupling during radial neuronal migration in the cerebral cortex and during glial migration (By similarity). Involved in telomere attachment to nuclear envelope in the prophase of meiosis implicating a SUN1/2:KASH5 LINC complex in which SUN1 and SUN2 seem to act at least partial



redundantly (By similarity). Required for gametogenesis and involved in selective gene expression of coding and non-coding RNAs needed for gametogenesis (By similarity). Helps to define the distribution of nuclear pore complexes (NPCs) (By similarity). Required for efficient localization of SYNE4 in the nuclear envelope (By similarity). May be involved in nuclear remodeling during sperm head formation in spermatogenesis (By similarity). May play a role in DNA repair by suppressing non-homologous end joining repair to facilitate the repair of DNA cross-links .

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images