



PDIR Rabbit pAb

Catalog No	YP-Ab-18457
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PDIA5 PDIR
Immunogen	Synthesized peptide derived from human PDIR
Specificity	This antibody detects endogenous levels of PDIR at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PDIA5;EC:5.3.4.1;Protein disulfide isomerase A5;Protein disulfide isomerase-related protein;Protein disulfide-isomerase A5;PDIA5 PDIR;Protein disulfide-isomerase A5 (Protein disulfide isomerase-related protein)
Observed Band	60 kDa
Calculated Molecular Weight	60 kDa
Cell Pathway	Endoplasmic reticulum lumen .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
