



UTP6 Rabbit pAb

Catalog No	YP-Ab-18459
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB
Gene Name	UTP6 C17orf40 HCA66 MHAT
Immunogen	Synthesized peptide derived from human UTP6
Specificity	This antibody detects endogenous levels of UTP6 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HCA66; C17orf40; UTP6; UTP6 C17orf40 HCA66 MHAT; U3 small nucleolar RNA-associated protein 6 homolog (Hepatocellular carcinoma-associated antigen 66) (Multiple hat domains protein)
Observed Band	70kDa
Calculated Molecular Weight	70kDa
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	Involved in nucleolar processing of pre-18S ribosomal RNA.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
