



SLITRK6 Mouse mAb

Catalog No	YP-Ab-18475
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB
Gene Name	SLITRK6
Immunogen	Synthesized peptide derived from human SLITRK6
Specificity	This antibody detects endogenous levels of SLITRK6 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DFNMYP; SLITRK6; SLIT and NTRK-like protein 6
Observed Band	120kDa
Calculated Molecular Weight	95kDa
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	In adult brain, highly expressed in putamen with no expression in cerebral cortex. Expressed in adult and fetal lung and fetal liver. Also expressed at high levels in some brain tumors including medulloblastomas and primitive neuroectodermal tumors.
Function	Regulator of neurite outgrowth required for normal hearing and vision.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images