



RBCK1 Rabbit pAb

Catalog No	YP-Ab-18496
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RBCK1 C20orf18 RNF54 UBCE7IP3 XAP3 XAP4
Immunogen	Synthesized peptide derived from human RBCK1
Specificity	This antibody detects endogenous levels of RBCK1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	XAP3;XAP4;HOIL1;PBMEI;PGBM1;RBCK2;RNF54;HOIL-1;ZRANB4;C20orf18;UBCE7IP3;RBCK1;RBCK1 C20orf18 RNF54 UBCE7IP3 XAP3 XAP4;RanBP-type and C3HC4-type zinc finger-containing protein 1 (HBV-associated factor 4) (Heme-oxidized IRP2 ubiquitin ligase 1) (HOIL-1) (Hepatitis B virus X-associated protein 4) (RING finger protein 54) (Ubiquitin-conjugating enzyme 7-interacting protein 3)
Observed Band	58kDa
Calculated Molecular Weight	58kDa
Cell Pathway	
Tissue Specificity	
Function	E3 ubiquitin-protein ligase, which accepts ubiquitin from specific E2 ubiquitin-conjugating enzymes, such as UBE2L3/UBCM4, and then transfers it to substrates . Functions as an E3 ligase for oxidized IREB2 and both heme and oxygen are necessary for IREB2 ubiquitination . Promotes ubiquitination of TAB2 and IRF3 and their degradation by the proteasome . Component of the LUBAC complex which conjugates linear ('Met-1'-linked) polyubiquitin chains to substrates and plays a key role in NF-kappa-B activation and regulation of inflammation . LUBAC conjugates linear polyubiquitin to IKBKG and RIPK1 and



is involved in activation of the canonical NF-kappa-B and the JNK signaling pathways . Linear ubiquitination mediated by the LUBAC complex interferes with TNF-induced cell death and thereby prevents inflammation . LUBAC is recruited to the TNF-R1 signaling complex (TNF-RSC) following polyubiquitination of TNF-RSC components by BIRC2 and/or BIRC3 and to conjugate linear polyubiquitin to IKBKG and possibly other components contributing to the stability of the complex . The LUBAC complex is also involved in innate immunity by conjugating linear polyubiquitin chains at the surface of bacteria invading the cytosol to form the ubiquitin coat surrounding bacteria . LUBAC is not able to initiate formation of the bacterial ubiquitin coat, and can only promote formation of linear polyubiquitins on pre-existing ubiquitin . The bacterial ubiquitin coat acts as an 'eat-me' signal for xenophagy and promotes NF-kappa-B activation . Together with OTULIN, the LUBAC complex regulates the canonical Wnt signaling during angiogenesis . Binds polyubiquitin of different linkage types .

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images