



ADH6 Rabbit pAb

Catalog No	YP-Ab-18521
Isotype	IgG
Reactivity	Human,Rat
Applications	WB
Gene Name	ADH6
Immunogen	Synthesized peptide derived from human ADH6
Specificity	This antibody detects endogenous levels of ADH6 at Human,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADH 5;ADH 6;Adh6;ADH6_HUMAN;Alcohol dehydrogenase 6 (class V);Alcohol dehydrogenase 6;Aldehyde reductase;Class V alcohol dehydrogenase 6;UDP GLUCURONOSYLTRANSFERASE 2B1 PRECURSOR;MICROSOMAL;Alcohol dehydrogenase 6 (EC 1.1.1.1)
Observed Band	40kD
Calculated Molecular Weight	39 kDa; 39kD(Calculated).
Cell Pathway	Cytoplasm.
Tissue Specificity	Stomach and liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
