



NGEF Rabbit pAb

Catalog No	YP-Ab-18530
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	NGEF
Immunogen	Synthesized peptide derived from human NGEF
Specificity	This antibody detects endogenous levels of NGEF at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EPHEXIN;ARHGEF27;NGEF;Ephexin-1 (Eph-interacting exchange protein) (Neuronal guanine nucleotide exchange factor)
Observed Band	82kDa
Calculated Molecular Weight	82kDa
Cell Pathway	Cytoplasm . Membrane . Cell projection, growth cone . Associated with membranes. Localizes to axonal growth cones (By similarity). .
Tissue Specificity	Highly expressed in brain specifically in caudate nucleus and to a lower extent in amygdala and hippocampus. Also detected in lung.
Function	Acts as a guanine nucleotide exchange factor (GEF) which differentially activates the GTPases RHOA, RAC1 and CDC42. Plays a role in axon guidance regulating ephrin-induced growth cone collapse and dendritic spine morphogenesis. Upon activation by ephrin through EPHA4, the GEF activity switches toward RHOA resulting in its activation. Activated RHOA promotes cone retraction at the expense of RAC1- and CDC42-stimulated growth cone extension (By similarity).
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images