



DEK Rabbit pAb

Catalog No	YP-Ab-18546
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	DEK
Immunogen	Synthesized peptide derived from human DEK
Specificity	This antibody detects endogenous levels of DEK at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D6S231E;DEK;Protein DEK
Observed Band	50kDa
Calculated Molecular Weight	43kDa
Cell Pathway	Nucleus . Enriched in regions where chromatin is decondensed or sparse in the interphase nuclei.
Tissue Specificity	Ubiquitous. Expressed at relatively high levels.
Function	Involved in chromatin organization.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
