



DLEU1 Rabbit pAb

Catalog No	YP-Ab-18551
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	DLEU1 LEU1 XTP6
Immunogen	Synthesized peptide derived from human DLEU1
Specificity	This antibody detects endogenous levels of DLEU1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DLEU1 LEU1 XTP6;Leukemia-associated protein 1 (Deleted in lymphocytic leukemia 1) (HBV X-transactivated gene 6 protein) (HBV XAg-transactivated protein 6)
Observed Band	9kD
Calculated Molecular Weight	3 kDa; 9kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	May act as a tumor suppressor.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
