



HYAL1 Rabbit pAb

Catalog No	YP-Ab-18569
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HYAL1 LUCA1
Immunogen	Synthesized peptide derived from human HYAL1
Specificity	This antibody detects endogenous levels of HYAL1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MPS9;NAT6;LUCA1;HYAL-1;HYAL1;HYAL1 LUCA1;Hyaluronidase-1 (Hyal-1) (Hyaluronoglucosaminidase-1) (Lung carcinoma protein 1) (LuCa-1)
Observed Band	48kDa
Calculated Molecular Weight	48kDa
Cell Pathway	
Tissue Specificity	
Function	May have a role in promoting tumor progression. May block the TGFB1-enhanced cell growth.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
