



FAM107A Rabbit pAb

Catalog No	YP-Ab-18604
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	FAM107A DRR1 TU3A
Immunogen	Synthesized peptide derived from human FAM107A
Specificity	This antibody detects endogenous levels of FAM107A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DRR1;TU3A;FAM107A;FAM107A DRR1 TU3A;Protein FAM107A (Down-regulated in renal cell carcinoma 1) (Protein TU3A)
Observed Band	18kDa
Calculated Molecular Weight	17kDa
Cell Pathway	Nucleus . Cytoplasm, cytoskeleton, stress fiber . Cell junction, focal adhesion . Cell projection, ruffle membrane . Cell junction, synapse . Colocalizes with F-actin and COMMD1 in the nucleus (PubMed:28604741). Colocalizes with MAP1A along actin stress fibers and membrane ruffles (PubMed:20543869). .
Tissue Specificity	Widely expressed (PubMed:10564580). Expressed in neurons (PubMed:20543869). Expressed in malignant glial tumors (PubMed:20543869). Expression is reduced or absent in a number of cancer cell lines (PubMed:10564580).
Function	Stress-inducible actin-binding protein that plays a role in synaptic and cognitive functions by modulating actin filamentous (F-actin) dynamics. Mediates polymerization of globular actin to F-actin. Also binds to, stabilizes and bundles F-actin. Involved in synaptic function by regulating neurite outgrowth in an actin-dependent manner and for the acquisition of hippocampus-dependent cognitive function, such as learning and long-term memory (By similarity). Plays a role in



the actin and microtubule cytoskeleton organization; negatively regulates focal adhesion (FA) assembly promoting malignant glial cell migration in an actin-, microtubule- and MAP1A-dependent manner . Also involved in neuroblastoma G1/S phase cell cycle progression and cell proliferation inhibition by stimulating ubiquitination of NF-kappa-B subunit RELA and NF-kappa-B degradation in a COMMD1- and actin-dependent manner . May play a role in tumor development .

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images