



PRM2 Rabbit pAb

Catalog No	YP-Ab-18631
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PRM2
Immunogen	Synthesized peptide derived from human PRM2
Specificity	This antibody detects endogenous levels of PRM2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000; IHC 1:100-1:300; ELISA 1:40000; IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Sperm histone HP2;Protamine-2;Basic nuclear protein HPS2;Basic nuclear protein HPS1;Basic nuclear protein HPI2;PRM2;Protamine-2 (Sperm histone P2) (Sperm protamine P2) [Cleaved into: Basic nuclear protein HPI1;Sperm histone HP4 (Sperm protamine P4);Sperm histone HP2 (Sperm protamine P2) (P2'); Sperm histone HP3 (P2'') (Sperm protamine P3)]
Observed Band	10-25 kDa
Calculated Molecular Weight	13 kDa
Cell Pathway	Nucleus . Chromosome .
Tissue Specificity	Testis.
Function	Protamines substitute for histones in the chromatin of sperm during the haploid phase of spermatogenesis. They compact sperm DNA into a highly condensed, stable and inactive complex.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images