



SCAMP3 Rabbit pAb

Catalog No	YP-Ab-18666
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB
Gene Name	SCAMP3 C1orf3 PROPIN1
Immunogen	Synthesized peptide derived from human SCAMP3
Specificity	This antibody detects endogenous levels of SCAMP3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1orf3; PROPIN1; Secretory carrier-associated membrane protein 3; SCAMP3 C1orf3 PROPIN1; Secretory carrier-associated membrane protein 3 (Secretory carrier membrane protein 3)
Observed Band	33 kDa
Calculated Molecular Weight	38 kDa
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Widely expressed, with highest expression in heart and skeletal muscle.
Function	Functions in post-Golgi recycling pathways. Acts as a recycling carrier to the cell surface.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
