



NPR3 Rabbit pAb

Catalog No	YP-Ab-18683
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	NPR3 ANPRC C5orf23 NPRC
Immunogen	Synthesized peptide derived from human NPR3
Specificity	This antibody detects endogenous levels of NPR3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NPRC;ANP-C;ANPRC;BOMOS;NPR-C;ANPR-C;GUCY2B;C5orf23;NPR3;NPR3 ANPRC C5orf23 NPRC;Atrial natriuretic peptide receptor 3 (Atrial natriuretic peptide clearance receptor) (Atrial natriuretic peptide receptor type C) (ANP-C) (ANPR-C) (NPR-C)
Observed Band	60kDa
Calculated Molecular Weight	60kDa
Cell Pathway	Membrane; Single-pass type I membrane protein.
Tissue Specificity	
Function	Receptor for the natriuretic peptide hormones, binding with similar affinities atrial natriuretic peptide NPPA/ANP, brain natriuretic peptide NPPB/BNP, and C-type natriuretic peptide NPPC/CNP. May function as a clearance receptor for NPPA, NPPB and NPPC, regulating their local concentrations and effects. May regulate diuresis, blood pressure and skeletal development. Does not have guanylate cyclase activity.
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images