



BEND5 Mouse pAb

Catalog No	YP-Ab-18690
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	BEND5 C1orf165
Immunogen	Synthesized peptide derived from human BEND5
Specificity	This antibody detects endogenous levels of BEND5 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BEN domain containing 5;BEND5;C1orf165;BEND5 C1orf165;BEN domain-containing protein 5
Observed Band	48 kDa
Calculated Molecular Weight	421 aa, 48 kDa
Cell Pathway	
Tissue Specificity	
Function	Acts as a transcriptional repressor .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
