



FUT8 Rabbit pAb

Catalog No	YP-Ab-18692
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	FUT8
Immunogen	Synthesized peptide derived from human FUT8
Specificity	This antibody detects endogenous levels of FUT8 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CDGF;CDGF1;FUT8;Alpha-(1;6)-fucosyltransferase (alpha1-6FucT) (Fucosyltransferase 8) (GDP-L-Fuc:N-acetyl-beta-D-glucosaminide alpha1;6-fucosyltransferase) (GDP-fucose--glycoprotein fucosyltransferase) (Glycoprotein 6-alpha-L-fucosyltransferase)
Observed Band	75kDa
Calculated Molecular Weight	67kDa
Cell Pathway	Golgi apparatus, Golgi stack membrane ; Single-pass type II membrane protein . Membrane-bound form in trans cisternae of Golgi. .
Tissue Specificity	
Function	Catalyzes the addition of fucose in alpha 1-6 linkage to the first GlcNAc residue, next to the peptide chains in N-glycans.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images