



# C2orf79/PTRHD1 Rabbit pA

|                                    |                                                                                                                                                                                                                                           |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-18709                                                                                                                                                                                                                               |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                       |
| <b>Reactivity</b>                  | Human, Mouse                                                                                                                                                                                                                              |
| <b>Applications</b>                | WB                                                                                                                                                                                                                                        |
| <b>Gene Name</b>                   | PTRHD1 C2orf79                                                                                                                                                                                                                            |
| <b>Immunogen</b>                   | Synthesized peptide derived from human C2orf79/PTRHD1                                                                                                                                                                                     |
| <b>Specificity</b>                 | This antibody detects endogenous levels of C2orf79/PTRHD1 at Human, Mouse                                                                                                                                                                 |
| <b>Formulation</b>                 |                                                                                                                                                                                                                                           |
| <b>Source</b>                      |                                                                                                                                                                                                                                           |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                     |
| <b>Dilution</b>                    | WB 1:500-2000                                                                                                                                                                                                                             |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                   |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                      |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                              |
| <b>Synonyms</b>                    | PTRHD1; C2orf79; Putative peptidyl-tRNA hydrolase PTRHD1; Peptidyl-tRNA hydrolase domain-containing protein 1; EC:3.1.1.29; PTRHD1 C2orf79; Putative peptidyl-tRNA hydrolase PTRHD1 (Peptidyl-tRNA hydrolase domain-containing protein 1) |
| <b>Observed Band</b>               | 140 aa, 16 kDa                                                                                                                                                                                                                            |
| <b>Calculated Molecular Weight</b> | 140 aa, 16 kDa                                                                                                                                                                                                                            |
| <b>Cell Pathway</b>                |                                                                                                                                                                                                                                           |
| <b>Tissue Specificity</b>          |                                                                                                                                                                                                                                           |
| <b>Function</b>                    |                                                                                                                                                                                                                                           |
| <b>Background</b>                  |                                                                                                                                                                                                                                           |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                                                                                                                      |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                                                         |



## Products Images

|  |
|--|
|  |
|  |
|  |
|  |
|  |
|  |
|  |