



TNS4 Rabbit pAb

Catalog No	YP-Ab-18772
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TNS4 CTEN PP14434
Immunogen	Synthesized peptide derived from human TNS4
Specificity	This antibody detects endogenous levels of TNS4 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C terminal tensin like;C terminal tensin like protein;C-terminal tensin-like protein; CTEN;FLJ1495;PP14434;TENS4_HUMAN;Tensin 4;Tensin-4;Tns4;TNS4 CTEN PP14434;Tensin-4 (C-terminal tensin-like protein)
Observed Band	79kD
Calculated Molecular Weight	77kD(Calculated).
Cell Pathway	Prostate and placenta. Down regulated in prostate cancer.
Tissue Specificity	
Function	Cell junction, focal adhesion . Cytoplasm, cytoskeleton .
Background	May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton (By similarity). May promote apoptosis, via its cleavage by caspase-3.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

