



PADI3/PAD3 Rabbit pAb

Catalog No	YP-Ab-18778
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PADI3 PDI3
Immunogen	Synthesized peptide derived from human PADI3/PAD3
Specificity	This antibody detects endogenous levels of PADI3/PAD3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.5.3.15;PAD3;PDI3;Peptidylarginine deiminase III;Protein-arginine deiminase type III;PADI3 PDI3;Protein-arginine deiminase type-3 (Peptidylarginine deiminase III) (Protein-arginine deiminase type III)
Observed Band	70-75 kDa
Calculated Molecular Weight	75kDa,664aa
Cell Pathway	Cytoplasm .
Tissue Specificity	Hair follicles, and epidermis at very low levels.
Function	Catalyzes the deimination of arginine residues of proteins.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

