



SCCPDH Rabbit pAb

Catalog No	YP-Ab-18810
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SCCPDH CGI-49
Immunogen	Synthesized peptide derived from human SCCPDH
Specificity	This antibody detects endogenous levels of SCCPDH at Human, Mouse,Rat
Formulation	
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NET11;CGI-49;SCCPDH;SCCPDH CGI-49;Saccharopine dehydrogenase-like oxidoreductase (EC 1.-.-.-)
Observed Band	45kDa
Calculated Molecular Weight	47kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

