



NKX3-2 Rabbit pAb

Catalog No	YP-Ab-18825
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	NKX3-2 BAPX1 NKX3B
Immunogen	Synthesized peptide derived from human NKX3-2
Specificity	This antibody detects endogenous levels of NKX3-2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SMMD;BAPX1;NKX3B;NKX3.2;NKX3-2;NKX3-2 BAPX1 NKX3B;Homeobox protein Nkx-3.2 (Bagpipe homeobox protein homolog 1) (Homeobox protein NK-3 homolog B)
Observed Band	40kDa/46kDa
Calculated Molecular Weight	35kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed at highest levels in cartilage, bone (osteosarcoma) and gut (small intestine and colon), whereas moderate expression is seen in trachea and brain. Expressed in visceral mesoderm and embryonic skeleton.
Function	Transcriptional repressor that acts as a negative regulator of chondrocyte maturation. PLays a role in distal stomach development; required for proper antral-pyloric morphogenesis and development of antral-type epithelium. In concert with GSC, defines the structural components of the middle ear; required for tympanic ring and gonium development and in the regulation of the width of the malleus (By similarity).
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images