



iASPP Rabbit pAb

Catalog No	YP-Ab-18848
Isotype	IgG
Reactivity	Human; Mouse;
Applications	WB
Gene Name	PPP1R13L IASPP NKIP1 PPP1R13BL RAI
Immunogen	Synthesized peptide derived from human iASPP
Specificity	This antibody detects endogenous levels of iASPP at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IASPP;inhibitor of apoptosis stimulating protein of p53;Inhibitor of ASPP protein; NFkB interacting protein 1;NKIP1;PPP1R13B-like protein;PPP1R13L;Protein iASPP;protein phosphatase 1 regulatory (inhibitor) subunit 13 like;RAI;RelA-associated inhibitor;PPP1R13L IASPP NKIP1 PPP1R13BL RAI;RelA-associated inhibitor (Inhibitor of ASPP protein) (Protein iASPP) (NFkB-interacting protein 1) (PPP1R13B-like protein)
Observed Band	
Calculated Molecular Weight	89kDa; 89kD(Calculated).
Cell Pathway	Cytoplasm . Nucleus . Predominantly cytoplasmic but also nuclear. .
Tissue Specificity	
Function	Regulator that plays a central role in regulation of apoptosis and transcription via its interaction with NF-kappa-B and p53/TP53 proteins. Blocks transcription of HIV-1 virus by inhibiting the action of both NF-kappa-B and SP1. Also inhibits p53/TP53 function, possibly by preventing the association between p53/TP53 and ASPP1 or ASPP2, and therefore suppressing the subsequent activation of apoptosis .



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images