



TRIM67 Rabbit pAb

Catalog No	YP-Ab-18849
Isotype	IgG
Reactivity	Human; Mouse;
Applications	WB
Gene Name	TRIM67 TNL
Immunogen	Synthesized peptide derived from human TRIM67
Specificity	This antibody detects endogenous levels of TRIM67 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TNL;TRI67_HUMAN;Trim67;TRIM9 like protein;TRIM9-like protein;Tripartite motif containing 67;Tripartite motif-containing protein 67;TRIM67 TNL;Tripartite motif-containing protein 67 (TRIM9-like protein)
Observed Band	
Calculated Molecular Weight	84kD(Calculated).
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton . Microtubule-associated. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images