



DTX3L Rabbit pAb

Catalog No	YP-Ab-18864
Isotype	IgG
Reactivity	Human; Mouse;
Applications	WB
Gene Name	DTX3L BBAP
Immunogen	Synthesized peptide derived from human DTX3L
Specificity	This antibody detects endogenous levels of DTX3L at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BBAP;RNF143;DTX3L;DTX3L BBAP;E3 ubiquitin-protein ligase DTX3L (B-lymphoma- and BAL-associated protein) (Protein deltex-3-like) (Rhysin-2) (Rhysin2)
Observed Band	84kDa
Calculated Molecular Weight	84kDa
Cell Pathway	Cytoplasm . Nucleus . Early endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Lysosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Translocates to the nucleus in response to IFNG or IFNB1 stimulation (PubMed:26479788). Localizes at sites of DNA damage in a PARP1-dependent manner (PubMed:23230272). Localization to early endosomes is increased upon CXCL12 stimulation where it co-localizes with ITCH, CXCL4, HGS and STAM (PubMed:24790097). A minor proportion localizes to lysosomes (PubMed:24790097). .
Tissue Specificity	
Function	E3 ubiquitin-protein ligase which, in association with ADP-ribosyltransferase PARP9, plays a role in DNA damage repair and in interferon-mediated antiviral responses . Monoubiquitinates several histones, including histone H2A, H2B, H3 and H4 . In response to DNA damage, mediates monoubiquitination of 'Lys-91'



of histone H4 (H4K91ub1) . The exact role of H4K91ub1 in DNA damage response is still unclear but it may function as a licensing signal for additional histone H4 post-translational modifications such as H4 'Lys-20' methylation (H4K20me) . PARP1-dependent PARP9-DTX3L-mediated ubiquitination promotes the rapid and specific recruitment of 53BP1/TP53BP1, UIMC1/RAP80, and BRCA1 to DNA damage sites . By monoubiquitinating histone H2B H2BC9/ H2BJ and thereby promoting chromatin remodeling, positively regulates STAT1-dependent interferon-stimulated gene transcription and thus STAT1-mediated control of viral replication . Independently of its catalytic activity, promotes the sorting of chemokine receptor CXCR4 from early endosome to lysosome following CXCL12 stimulation by reducing E3 ligase ITCH activity and thus ITCH-mediated ubiquitination of endosomal sorting complex required for transport ESCRT-0 components HGS and STAM . In addition, required for the recruitment of HGS and STAM to early endosomes . In association with PARP9, plays a role in antiviral responses by mediating 'Lys-48'-linked ubiquitination of encephalomyocarditis virus (EMCV) and human rhinovirus (HRV) C3 proteases and thus promoting their proteosomal-mediated degradation .

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images