



NXT2 Rabbit pAb

Catalog No	YP-Ab-18890
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	NXT2 BM-025 DC9
Immunogen	Synthesized peptide derived from human NXT2
Specificity	This antibody detects endogenous levels of NXT2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P15-2;NXT2;NXT2 BM-025 DC9;NTF2-related export protein 2 (Protein p15-2)
Observed Band	16kDa
Calculated Molecular Weight	16kDa
Cell Pathway	Nucleus. Cytoplasm. Shuttles between the nucleus and the cytoplasm.
Tissue Specificity	
Function	Regulator of protein export for NES-containing proteins. Also plays a role in mRNA nuclear export.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images