



SLMAP Rabbit pAb

Catalog No	YP-Ab-18920
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SLMAP KIAA1601 SLAP UNQ1847/PRO3577
Immunogen	Synthesized peptide derived from human SLMAP
Specificity	This antibody detects endogenous levels of SLMAP at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLAP;SLMAP;SLMAP KIAA1601 SLAP UNQ1847/PRO3577;Sarcolemmal membrane-associated protein (Sarcolemmal-associated protein)
Observed Band	110kDa
Calculated Molecular Weight	95kDa
Cell Pathway	Cell membrane, sarcolemma ; Single-pass type IV membrane protein . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Membrane-associated. Distributed in the transverse tubules and near the junctional sarcoplasmic reticulum. Detected along the Z- and M-lines in cardiomyocytes. Centrosome. Localizes to the centrosomes in a microtubule-dependent manner (By similarity). .
Tissue Specificity	
Function	May play a role during myoblast fusion.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images