



GIP Rabbit pAb

Catalog No	YP-Ab-18936
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	GIP
Immunogen	Synthesized peptide derived from human GIP
Specificity	This antibody detects endogenous levels of GIP at Human, Mouse,Rat
Formulation	
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gastric Inhibitory Peptide;Gastric inhibitory polypeptide;Gastric inhibitory polypeptide precursor;GIP;GIP_HUMAN;Glucose dependent insulinotropic polypeptide;Glucose-dependent insulinotropic polypeptide;Incretin hormone; Gastric inhibitory polypeptide (GIP) (Glucose-dependent insulinotropic polypeptide) (Incretin hormone)
Observed Band	
Calculated Molecular Weight	17kD(Calculated).
Cell Pathway	Secreted.
Tissue Specificity	
Function	Potent stimulator of insulin secretion and relatively poor inhibitor of gastric acid secretion.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images