



C13orf1 Rabbit pAb

Catalog No	YP-Ab-18938
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SPRYD7 C13orf1 CLLD6
Immunogen	Synthesized peptide derived from human C13orf1
Specificity	This antibody detects endogenous levels of C13orf1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPRYD7;SPRY7;SPRY domain-containing protein 7;CLLD6;SPRYD7 C13orf1 CLLD6;SPRY domain-containing protein 7 (Chronic lymphocytic leukemia deletion region gene 6 protein) (CLL deletion region gene 6 protein)
Observed Band	17 kDa, 2 kDa
Calculated Molecular Weight	22kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images