



PLA2G12A Rabbit pAb

Catalog No	YP-Ab-18942
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	PLA2G12A PLA2G12 FKSG38 UNQ2519/PRO6012
Immunogen	Synthesized peptide derived from human PLA2G12A
Specificity	This antibody detects endogenous levels of PLA2G12A at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC 3.1.1.4;Group XII secreted phospholipase A2;Group XIIA secreted phospholipase A2;Group XIIA secretory phospholipase A2;Group XIIA secretory phospholipase A2 precursor;GXII;GXII sPLA2;PG12A_HUMAN;phosphatidylcholine 2 acylhydrolase 12A;Phosphatidylcholine 2 acylhydrolase GXII;Phosphatidylcholine 2-acylhydrolase 12A;Phospholipase A2 group XIIA;PLA2G12;Pla2g12a;ROSSY;sPLA2 XII;sPLA2-XII;PLA2G12A PLA2G12 FKSG38 UNQ2519/PRO6012;Group XIIA secretory phospholipase A2 (GXII sPLA2) (sPLA2-XII) (Phosphatidylcholine 2-acylhydrolase 12A)
Observed Band	
Calculated Molecular Weight	21 kDa; 21kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	PA2 catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides. Does not exhibit detectable activity toward sn-2-arachidonoyl- or linoleoyl-phosphatidylcholine or -phosphatidylethanolamine.
Background	



**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images