



WSB1 Rabbit pAb

Catalog No	YP-Ab-18949
Isotype	IgG
Reactivity	Human,Mouse
Applications	
Gene Name	
Immunogen	Synthesized peptide derived from human WSB1
Specificity	This antibody detects endogenous levels of WSB1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SWIP1;WSB-1;WSB1;WSB1 SWIP1;WD repeat and SOCS box-containing protein 1 (WSB-1) (SOCS box-containing WD protein SWIP-1)
Observed Band	47kDa
Calculated Molecular Weight	47kDa
Cell Pathway	
Tissue Specificity	
Function	Probable substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Recognizes type II iodothyronine deiodinase/DIO2. Confers constitutive instability to HIPK2 through proteasomal degradation.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images