



Cyclon/CCDC86 Rabbit pAb

Catalog No	YP-Ab-18970
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	CCDC86 CYCLON
Immunogen	Synthesized peptide derived from human Cyclon/CCDC86
Specificity	This antibody detects endogenous levels of Cyclon/CCDC86 at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CYCLON;CCDC86;Cytokine-induced protein with coiled-coil domain;Coiled-coil domain-containing protein 86;CCDC86 CYCLON;Coiled-coil domain-containing protein 86 (Cytokine-induced protein with coiled-coil domain)
Observed Band	65 kDa
Calculated Molecular Weight	40 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images