



ADORA2B Rabbit pAb

Catalog No	YP-Ab-18977
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	ADORA2B
Immunogen	Synthesized peptide derived from human ADORA2B
Specificity	This antibody detects endogenous levels of ADORA2B at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADORA2;ADORA2B;Adenosine receptor A2b
Observed Band	36kDa/30-50 kDa
Calculated Molecular Weight	36kDa
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	Receptor for adenosine. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images