



TSPAN2 Rabbit pAb

Catalog No	YP-Ab-19022
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TSPAN2
Immunogen	Synthesized peptide derived from human TSPAN2
Specificity	This antibody detects endogenous levels of TSPAN2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	6330415F13Rik;NET3;Tetraspan NET 3;tetraspanin 2;TSN2;TSPAN 2;TSPAN2; Tetraspanin-2 (Tspan-2) (Tetraspan NET-3)
Observed Band	30-35 kDa
Calculated Molecular Weight	221 aa, 24 kDa
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	May play a role in signalling in oligodendrocytes in the early stages of their terminal differentiation into myelin-forming glia and may also function in stabilizing the mature sheath.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images