



SKAP2 Rabbit pAb

Catalog No	YP-Ab-19042
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SKAP2 PRAP RA70 SAPS SCAP2 SKAP55R
Immunogen	Synthesized peptide derived from human SKAP2
Specificity	This antibody detects endogenous levels of SKAP2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRAP;RA70;SAPS;SCAP2;SKAP55R;SKAP-HOM;SKAP2;SKAP2 PRAP RA70 SAPS SCAP2 SKAP55R;Src kinase-associated phosphoprotein 2 (Pyk2/RAFTK-associated protein) (Retinoic acid-induced protein 70) (SKAP55 homolog) (SKAP-55HOM) (SKAP-HOM) (Src family-associated phosphoprotein 2) (Src kinase-associated phosphoprotein 55-related protein) (Src-associated adapter protein with PH and SH3 domains)
Observed Band	50kDa
Calculated Molecular Weight	41kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitously expressed. Present in platelets (at protein level).
Function	May be involved in B-cell and macrophage adhesion processes. In B-cells, may act by coupling the B-cell receptor (BCR) to integrin activation. May play a role in src signaling pathway.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images