



HSDL2 Rabbit pAb

Catalog No	YP-Ab-19049
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HSDL2 C9orf99
Immunogen	Synthesized peptide derived from human HSDL2
Specificity	This antibody detects endogenous levels of HSDL2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C9orf99;EC:1.-.-.-;Hydroxysteroid dehydrogenase-like protein 2;SDR13C1;Short chain dehydrogenase/reductase family 13C member 1;HSDL2 C9orf99; Hydroxysteroid dehydrogenase-like protein 2 (EC 1.-.-.-)
Observed Band	45 kDa
Calculated Molecular Weight	45 kDa
Cell Pathway	Peroxisome .
Tissue Specificity	Ubiquitous.
Function	Has apparently no steroid dehydrogenase activity.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images