



PRPSAP1 Rabbit pAb

Catalog No	YP-Ab-19065
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PRPSAP1
Immunogen	Synthesized peptide derived from human PRPSAP1
Specificity	This antibody detects endogenous levels of PRPSAP1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	39 kDa phosphoribosypyrophosphate synthase-associated protein;PAP39; Phosphoribosyl pyrophosphate synthase-associated protein 1;PRPP synthase-associated protein 1;PRPSAP1;Phosphoribosyl pyrophosphate synthase-associated protein 1 (PRPP synthase-associated protein 1) (39 kDa phosphoribosypyrophosphate synthase-associated protein) (PAP39)
Observed Band	39 kDa
Calculated Molecular Weight	356 aa, 39 kDa
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	Seems to play a negative regulatory role in 5-phosphoribose 1-diphosphate synthesis.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images