



HAGH Rabbit pAb

Catalog No	YP-Ab-19087
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HAGH GLO2 HAGH1
Immunogen	Synthesized peptide derived from human HAGH. AA range:1-30
Specificity	This antibody detects endogenous levels of HAGH at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GLO2;GLX2;GLXII;HAGH1;HAGH;HAGH GLO2 HAGH1; Hydroxyacylglutathione hydrolase;mitochondrial (Glyoxalase II) (Glx II)
Observed Band	34kDa
Calculated Molecular Weight	34kDa
Cell Pathway	[Isoform 1]: Mitochondrion matrix .; [Isoform 2]: Cytoplasm .
Tissue Specificity	Expressed in liver and kidney.
Function	Thiolesterase that catalyzes the hydrolysis of S-D-lactoyl-glutathione to form glutathione and D-lactic acid.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images