



FMOD Rabbit pAb

Catalog No	YP-Ab-19112
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	FMOD FM SLRR2E
Immunogen	Synthesized peptide derived from human Fibromodulin
Specificity	This antibody detects endogenous levels of Fibromodulin at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FMOD FM SLRR2E;Fibromodulin (FM) (Collagen-binding 59 kDa protein) (Keratan sulfate proteoglycan fibromodulin) (KSPG fibromodulin)
Observed Band	
Calculated Molecular Weight	41kD
Cell Pathway	Secreted, extracellular space, extracellular matrix.
Tissue Specificity	
Function	Affects the rate of fibrils formation. May have a primary role in collagen fibrillogenesis (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images