



RAB3C Rabbit pAb

Catalog No	YP-Ab-19129
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RAB3C
Immunogen	Synthesized peptide derived from human RAB3C
Specificity	This antibody detects endogenous levels of RAB3C at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RAB3C;RAB3C member RAS oncogene family;RAB3C_HUMAN;Ras related protein Rab3C;Ras-related protein Rab-3C
Observed Band	25kD
Calculated Molecular Weight	28 kDa; 26kD(Calculated).
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	Expressed in brain, placenta and lung.
Function	Protein transport. Probably involved in vesicular traffic (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images