



SLC25A20 Rabbit pAb

Catalog No	YP-Ab-19153
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SLC25A20 CAC CACT
Immunogen	Synthesized peptide derived from human SLC25A20
Specificity	This antibody detects endogenous levels of SLC25A20 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAC;CACT;SLC25A20;SLC25A20 CAC CACT;Mitochondrial carnitine/acylcarnitine carrier protein (Carnitine/acylcarnitine translocase) (CAC) (Solute carrier family 25 member 20)
Observed Band	33kDa
Calculated Molecular Weight	33kDa
Cell Pathway	Mitochondrion inner membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	Mediates the transport of acylcarnitines of different length across the mitochondrial inner membrane from the cytosol to the mitochondrial matrix for their oxidation by the mitochondrial fatty acid-oxidation pathway.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

