



TRIM26 Rabbit pAb

Catalog No	YP-Ab-19163
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TRIM26 RNF95 ZNF173
Immunogen	Synthesized peptide derived from human TRIM26
Specificity	This antibody detects endogenous levels of TRIM26 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Acid finger protein;AFP;RING finger protein 95;RNF95;TRIM26;tripartite motif containing 26;Zinc finger protein 173;ZNF173;TRIM26 RNF95 ZNF173;Tripartite motif-containing protein 26 (Acid finger protein) (AFP) (RING finger protein 95) (Zinc finger protein 173)
Observed Band	62 kDa
Calculated Molecular Weight	62 kDa
Cell Pathway	Cytoplasm . Nucleus . Viral infection mediates TRIM26 nuclear translocation. .
Tissue Specificity	
Function	E3 ubiquitin-protein ligase which regulates the IFN-beta production and antiviral response downstream of various DNA-encoded pattern-recognition receptors (PRRs). Promotes nuclear IRF3 ubiquitination and proteasomal degradation. Bridges together TBK1 and NEMO during the innate response to viral infection leading to the activation of TBK1.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images