



PYCR2 Rabbit pAb

Catalog No	YP-Ab-19172
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PYCR2
Immunogen	Synthesized peptide derived from human PYCR2
Specificity	This antibody detects endogenous levels of PYCR2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HLD10;P5CR2;PYCR2;Pyrroline-5-carboxylate reductase 2 (P5C reductase 2) (P5CR 2) (EC 1.5.1.2)
Observed Band	35kDa
Calculated Molecular Weight	34kDa
Cell Pathway	Cytoplasm . Mitochondrion .
Tissue Specificity	Detected in erythrocytes (at protein level) (PubMed:2722838, PubMed:6894153) . Expressed in fetal brain (PubMed:25865492).
Function	Housekeeping enzyme that catalyzes the last step in proline biosynthesis. In some cell types, such as erythrocytes, its primary function may be the generation of NADP(+). Can utilize both NAD and NADP. Has higher affinity for NADP, but higher catalytic efficiency with NADH . Involved in cellular response to oxidative stress .
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

