



RPL8 Rabbit pAb

Catalog No	YP-Ab-19179
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RPL8
Immunogen	Synthesized peptide derived from human RPL8
Specificity	This antibody detects endogenous levels of RPL8 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	L8;uL2;RPL8;60S ribosomal protein L8
Observed Band	28kDa
Calculated Molecular Weight	28kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	Component of the large ribosomal subunit.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

