



NPEPL1 Rabbit pAb

Catalog No	YP-Ab-19182
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	NPEPL1 KIAA1974
Immunogen	Synthesized peptide derived from human NPEPL1
Specificity	This antibody detects endogenous levels of NPEPL1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	aminopeptidase like 1;KIAA1974;NPEPL1;Probable aminopeptidase NPEPL1; NPEPL1 KIAA1974;Probable aminopeptidase NPEPL1 (Aminopeptidase-like 1)
Observed Band	52 kDa
Calculated Molecular Weight	523 aa, 56 kDa
Cell Pathway	Ubiquitously expressed.
Tissue Specificity	
Function	Probably catalyzes the removal of unsubstituted N-terminal amino acids from various peptides.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images