



SIAE Rabbit pAb

Catalog No	YP-Ab-19193
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SIAE YSG2
Immunogen	Synthesized peptide derived from human SIAE
Specificity	This antibody detects endogenous levels of SIAE at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LSE;AIS6;CSEC;YSG2;CSE-C;SIAE;SIAE YSG2;Sialate O-acetylerase (H-Lse) (Sialic acid-specific 9-O-acetylerase)
Observed Band	58kDa
Calculated Molecular Weight	58kDa
Cell Pathway	Lysosome . Secreted.
Tissue Specificity	Widely expressed with high expression in the testis, prostate, and colon.
Function	Catalyzes the removal of O-acetyl ester groups from position 9 of the parent sialic acid, N-acetylneuraminic acid.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

