



SAA3 Rabbit pAb

Catalog No	YP-Ab-19208
Isotype	IgG
Reactivity	Mouse
Applications	WB
Gene Name	SAA3
Immunogen	Synthesized peptide derived from Mouse SAA3 AA range: 30-110
Specificity	This antibody detects endogenous levels of Mouse SAA3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	I7R3;Saa-3;SAA3;Serum amyloid A-3 protein
Observed Band	14kDa
Calculated Molecular Weight	13kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

