



Hemopexin Rabbit pAb

Catalog No	YP-Ab-19227
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HPX
Immunogen	Synthesized peptide derived from human Hemopexin
Specificity	This antibody detects endogenous levels of Hemopexin at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Beta 1B glycoprotein;hemopexin;HPX;Hemopexin (Beta-1B-glycoprotein)
Observed Band	65-75 kDa
Calculated Molecular Weight	254 aa, 29 kDa
Cell Pathway	Secreted.
Tissue Specificity	Expressed by the liver and secreted in plasma.
Function	Binds heme and transports it to the liver for breakdown and iron recovery, after which the free hemopexin returns to the circulation.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images